

CAV1 Antibody



PACO13717

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

CAV1

Uniprot

Q03135

Synonyms:

caveolin 1, caveolae protein, 22kDa; CAV1; BSCL3; CAV; CGL3; MSTP085; VIP21 ;

Immunogen:

Human CAV1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A