

RPL37A Antibody



PACO13731

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

RPL37A

Uniprot

P61513

Synonyms:

ribosomal protein L37a; RPL37A; MGC74786 ;

Immunogen:

Human RPL37A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A