

TNS4 Antibody



PACO13745

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TNS4

Uniprot

Q8IZW8

Synonyms:

tensin 4; TNS4; CTEN; FLJ14950; PP14434 ;

Immunogen:

Human TNS4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A