

WEE2 Antibody



PACO13751

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

WEE2

Uniprot

P0C1S8

Synonyms:

WEE1 homolog 2 (S. pombe); WEE2; FLJ16107; WEE1B ;

Immunogen:

Human WEE2.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A