

WFS1 Antibody



PACO13752

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

WFS1

Uniprot

O76024

Synonyms:

Wolfram syndrome 1 (wolframin); WFS1; FLJ51211; WFRS; WFS; WOLFRAMIN ;

Immunogen:

Human WFS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A