

# PDE9A Antibody



PACO13758

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PDE9A

**Uniprot**

O76083

**Synonyms:**

PDE9A, High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A

**Immunogen:**

Human CD99.

**Storage:**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Product Images

---

N/A

N/A