

PACO17121

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:5000, IHC:1:25-1:100

Protein Background:

SNAI3 is a member of the SNAIL gene family, named for the Drosophila snail gene, which plays roles in mesodermal formation during embryogenesis.

Gene ID:

SNAI3

Uniprot

Q3KNW1

Synonyms:

snail family zinc finger 3


Immunogen:

Fusion protein of human SNAI3.

Storage:

-20° C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using PACO17121(SNAI3 Antibody) at dilution 1/25, on the right is treated with fusion protein. (Original magnification: x—200).