

# STXBP4 Antibody



PACO17195

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:5000, IHC:1:25-1:100

**Protein Background:**

**Gene ID:**

STXBP4

**Uniprot**

Q6ZWJ1

**Synonyms:**

syntaxin binding protein 4

**Immunogen:**


Fusion protein of human STXBP4.

**Storage:**

-20&deg; C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

## Product Images

---



The image on the left is immunohistochemistry of paraffin-embedded Human breast cancer tissue using PACO17195(STXBP4 Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: x—200).