

# THAP4 Antibody



PACO17230

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:5000, WB:1:500-1:2000,  
IHC:1:25-1:100

**Protein Background:**

THAP domain-containing protein 4 is a protein that in humans is encoded by the THAP4 gene.

**Gene ID:**

THAP4

**Uniprot**

Q8WY91

**Synonyms:**

THAP domain containing 4

**Immunogen:**


Fusion protein of human THAP4.

**Storage:**


-20°C; C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

## Product Images

---



Gel: 8%SDS-PAGE, Lysate: 80  $\mu$ g, Lane: 293T cell, Primary antibody: PACO17230(THAP4 Antibody) at dilution 1/250 dilution, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 10 seconds.



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using PACO17230(THAP4 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x—200).