

TMEM161A Antibody



PACO17273

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:5000, IHC:1:25-1:100

Protein Background:

May play a role in protection against oxidative stress. Overexpression leads to reduced levels of oxidant-induced DNA damage and apoptosis.

Gene ID:

TMEM161A

Uniprot

Q9NX61

Synonyms:

transmembrane protein 161A


Immunogen:

Fusion protein of human TMEM161A.

Storage:

-20° C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using PACO17273(TMEM161A Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x—200).