

TLR1 Antibody



PACO21065

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:**Gene ID:**

TLR1

Uniprot

Q15399

Synonyms:

TLR1; CD281; DKFZp547I0610; DKFZp564I0682; KIAA0012; MGC104956; MGC126311;
MGC126312; TIL; rsc786


Immunogen:

Recombinant protein of human TLR1.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of Spleen cell and Liver cell using TLR1 antibody.