

PDGFRA Antibody



PACO21249

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:**Gene ID:**

PDGFRA

Uniprot

P16234

Synonyms:

PDGFRA; CD140A; MGC74795; PDGFR2; Rhe-PDGFRA


Immunogen:

Recombinant protein of human PDGFRA.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of HepG2 cell lines, using PDGFRA antibody.