

TNFAIP3 Antibody



PACO21258

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:

Gene ID:

TNFAIP3

Uniprot

P21580

Synonyms:


TNFAIP3; A20; MGC104522; MGC138687; MGC138688; OTUD7C; TNFA1P2

Immunogen:

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of HeLa cell line, using TNFAIP3 antibody.