

# LETMD1 Antibody



PACO21269

---

## Product Information

**Size:**

100ul(100ug)

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:50-1:200

**Protein Background:**

This gene encodes a mitochondrial outer membrane protein. It has a potential role in tumorigenesis, which may result from negative regulation of the p53 tumor suppressor gene. Alternatively spliced transcript variants have been noted for this gene. [provided by RefSeq, Aug 2011].

**Gene ID:**

LETMD1

**Uniprot**

Q6P1Q0

**Synonyms:**

LETMD1; HCCR1; HCCR-1; HCCR-2; HCRR-2; 1110019O13Rik

**Immunogen:**


Recombinant protein of human LETMD1.

**Storage:**

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot analysis of extracts of PANC1 cell lines, using LETMD1 antibody.