

Scn4b Antibody



PACO21276

Product Information

Size:

100ul(100ug)

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:50-1:200

Protein Background:**Gene ID:**

Scn4b

Uniprot

Q7M729

Synonyms:

SCN4B; Gm1471


Immunogen:

Recombinant protein of mouse SCN4B.

Storage:

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of U-87 cell lines, using SCN4B antibody.