

# TGFB1 Antibody



PACO21329

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:50-1:200

**Protein Background:**

Multifunctional protein that controls proliferation, differentiation and other functions in many cell types. Many cells synthesize TGFB1 and have specific receptors for it. It positively and negatively regulates many other growth factors. .

**Gene ID:**

TGFB1

**Uniprot**

P01137

**Synonyms:**

CED; diaphyseal dysplasia 1, progressive; DPD1; TGF-beta 1; TGFB. TGF-BETA

**Immunogen:**


The antiserum was produced against Synthesized peptide derived from human TGF beta.

**Storage:**

1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol

## Product Images

---



Western blot analysis of extracts from HepG2 cells, using TGF beta Antibody. The lane on the left is treated with the synthesized peptide.