

BBC3 Antibody



PACO21344

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Gene ID:

BBC3

Uniprot

Q9BXH1

Synonyms:


BBC3, JFY1, JFY-1,

Immunogen:

Storage:

Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 and MCF-7 cells using Puma Antibody.