

# PCNA Antibody



PACO21345

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:2000,  
IHC:1:50-1:200

**Protein Background:**

**Gene ID:**

PCNA

**Uniprot**

P12004

**Synonyms:**

Proliferating cell nuclear antigen

**Immunogen:**


Fusion protein corresponding to a region derived from 9-252 amino acid, of Human Proliferating cell nuclear antigen.

**Storage:**


Purified Rabbit polyclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Product Images

---



The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using 31002 (PCNA Antibody) at dilution 1/25, on the right is treated with the fusion protein.



Gel: 10%SDS-PAGE Lysate: 40 ug Hela cell lysate Primary antibody: 1/500 dilution Secondary antibody: Donkey anti Rabbit IgG - H&L (HRP) at 1/3000 dilution Exposure time: 1 minute.