

VAV1 (Ab-174) Antibody



PACO21387

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Couples tyrosine kinase signals with the activation of the Rho/Rac GTPases, thus leading to cell differentiation and/or proliferation. Amarasinghe GK, et al.

Gene ID:

VAV1

Uniprot

P15498

Synonyms:

VAV


Immunogen:

Peptide sequence around aa. 172~176 (E-I-Y-E-D) derived from Human Vav.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from JK cells using Vav(Ab-174) Antibody.