

CTNNA1 (Ab-641) Antibody



PACO21421

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

CTNNA1

Uniprot

P35221

Synonyms:

Cadherin-associated protein; alpha E-catenin; NY-REN-13 antigen

Applications:

ELISA, WB

Immunogen:

Peptide sequence around aa. 639~643 (D-D-S-D-F) derived from Human a-catenin.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa and 293 cells using α -catenin(Ab-641) antibody.