

TFE3 Antibody



PACO21453

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

TFE3

Uniprot

P19532

Synonyms:

BHLHE33; TFEA; RCCP2


Immunogen:

Peptide sequence around aa.303~307 (L-D-V-Y-S) derived from Human TFE3.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells using TFE3 Antibody.