

GATA4 Antibody



PACO21482

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Transcriptional activator. Binds to the consensus sequence 5'-AGATAG-3'. Acts as a transcriptional activator of ANF in cooperation with NKX2-5 Hirayama-Yamada K., Kamisago M., Akimoto K., Aotsuka H., Nakamura Y., Tomita H., Furutani M., Matsuoka R. Am. J. Med. Genet. A 135:47-52(2005).

Gene ID:

GATA4

Uniprot

P43694

Synonyms:

MGC126629


Immunogen:

Peptide sequence around aa.434~438(D-S-H-G-D) derived from Human GATA-4.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extract from HeLa and 293 cells using GATA-4 Antibody.