

NACC1 Antibody



PACO21523

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

NACC1

Uniprot

Q96RE7

Synonyms:

BEND8; NACC1; NAC-1; BTBD14B; FLJ37383


Immunogen:

Peptide sequence around aa.519~523(A-G-P-S-A) derived from Human NAC1.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Rat and Mouse brain tissue using NAC1 Antibody.