

GFI1B Antibody



PACO21553

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:

Gene ID:

GFI1B

Uniprot

Q5VTD9

Synonyms:

Gfi-1b


Immunogen:

Peptide sequence around aa.322~326(H-R-E-S-Q) derived from Human GFI1b.

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Raji cells using GFI1b Antibody.