

HNF1B Antibody



PACO21556

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000

Protein Background:**Gene ID:**

HNF1B

Uniprot

P35680

Synonyms:

FJHN; HNF2; LFB3; TCF2; HPC11


Immunogen:

Peptide sequence around aa.252~256(R-Q-K-N-P) derived from Human HNF-1b(TCF-2).

Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa, K562 and 3T3 cells using HNF-1 β (TCF-2) Antibody.