

PPP1R2 (Ab-44) Antibody



PACO21630

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Inhibitor of protein-phosphatase 1. Helps N. R. , FEBS Lett. 340:93-98(1994).

Gene ID:

PPP1R2

Uniprot

P41236

Synonyms:

I-2; IPP-2; IPP2; Protein phosphatase inhibitor 2;


Immunogen:

Synthesized peptide derived from Internal of human PPP1R2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COLO cells, HepG2 cells and HUVEC cells, using PPP1R2 (Ab-44) antibody.