

COL18A1 Antibody



PACO21651

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

COL18A1

Uniprot

P39060

Synonyms:

COL18A1; ENDOSTATIN; collagen type XVIII alpha 1;

Applications:

ELISA, IHC

Immunogen:

Synthesized peptide derived from human endostatin.

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue, using Endostatin antibody.

