

FBN1 Antibody



PACO21652

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

FBN1

Uniprot

P35555

Synonyms:

MASS; MFS1; OCTD; SGS; WMS

Applications:

ELISA, IHC

Immunogen:

Synthesized peptide derived from human Fibrillin-1.

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue, using Fibrillin-1 antibody.

