

PACO21673

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Part of the TFIID complex, a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. Potentiates transcriptional activation by the AF-2S of the retinoic acid, vitamin D3 and thyroid hormone. Deloukas P. , Nature 414:865-871(2001).

Gene ID:

TAF4

Uniprot

O00268

Synonyms:

Transcription initiation factor TFIID subunit 4; TBP-associated factor 4; Transcription initiation factor TFIID 135 kDa subunit; TAF(II)135; TAFII-135


Immunogen:

Synthesized peptide derived from internal of human TAF4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from NIH-3T3 cells, using TAF4 antibody.