

PFDN1 Antibody



PACO21697

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:**Protein Background:****Gene ID:**

PFDN1

Uniprot

O60925

Synonyms:

Prefoldin subunit 1; PFDN1; PFD1;

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human PFDN1.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from RAW264.7 cells, using PFDN1 antibody.