

CIB2 Antibody



PACO21704

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

CIB2

Uniprot

O75838

Synonyms:

Calcium and integrin-binding family member 2; Kinase-interacting protein 2; KIP 2; CIB2; KIP2

Immunogen:

Synthesized peptide derived from internal of human CIB2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using CIB2 antibody.

