

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Belongs to the CDK2AP family. Kalnine N. , Submitted (MAY-2003) to the EMBL/GenBank/DDBJ databases.

Gene ID:

CDK2AP2

Uniprot

O75956

Synonyms:

Cyclin-dependent kinase 2-associated protein 2; CDK2-associated protein 2; DOC-1-related protein; DOC-1R; CDK2AP2

Immunogen:

Synthesized peptide derived from internal of human CDKA2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using CDKA2 antibody.

