

MAP3K6 Antibody



PACO21709

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Component of a protein kinase signal transduction cascade. Activates the JNK, but not ERK or p38 kinase pathways. Gregory S. G. , Nature 441:315-321(2006).

Gene ID:

MAP3K6

Uniprot

O95382

Synonyms:

EC 2.7.11.25; M3K6; MAPKKK6; Mitogen-activated protein kinase kinase kinase 6; kinase MAP3K6

Immunogen:

Synthesized peptide derived from internal of human MAP3K6.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using MAP3K6 antibody.

