

FGR Antibody



PACO21723

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

FGR

Uniprot

P09769

Synonyms:

C-FGR; EC 2.7.10.2; FGR protein; P55-FGR; Proto-oncogene tyrosine-protein kinase FGR

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human FGR.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells and COLO205 cells, using FGR antibody.