

POLA1 Antibody



PACO21724

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

POLA1

Uniprot

P09884

Synonyms:

DNA pol alpha; DNA polymerase alpha catalytic subunit; DPOA; DPOLA; EC 2.7.7.7

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from N-terminal of human DNA Polymerase α .


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, 293 cells and HepG2 cells, using DNA Polymerase alpha antibody.