

HOXD12 Antibody



PACO21775

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

HOXD12

Uniprot

P35452

Synonyms:

Homeobox protein Hox-D12; Hox-4H; HOXD12; HOX4H;


Immunogen:

Synthesized peptide derived from internal of human HOXD12.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562, JK and 3T3 cells using HOXD12 antibody.