

ZNF75D Antibody



PACO21794

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May be involved in transcriptional regulation. Villa A. , Genomics 13:1231-1236(1992).

Gene ID:

ZNF75D

Uniprot

P51815

Synonyms:

Zinc finger protein 75; Zinc finger protein 82; ZNF75; ZNF82;


Immunogen:

Synthesized peptide derived from internal of human ZNF75.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells, using ZNF75 antibody.