

# AFF1 Antibody



PACO21795

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

**Protein Background:**

**Gene ID:**

AFF1

**Uniprot**

P51825

**Synonyms:**

AF4; AFF1; FEL protein; MLLT2; Proto-oncogene AF4

**Applications:**

ELISA, IHC

**Immunogen:**

Synthesized peptide derived from N-terminal of human AF4.

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using AF4 antibody.

