

U69 Antibody



PACO21799

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Phosphorylates the antiviral nucleoside analog ganciclovir By similarity.

Gene ID:

U69

Uniprot

P52446

Synonyms:

Probable ganciclovir kinase; EC=2.7.1;

Immunogen:

Synthesized peptide derived from internal of human GCVK_HHV6Z.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

N/A

N/A