

TAF1A Antibody



PACO21858

Product Information

Size:

100ul

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

TAF1A

Uniprot

Q15573

Synonyms:

TATA box-binding protein-associated factor 1A; TBP-associated factor 1A; TBP-associated factor RNA polymerase I 48 kDa; TAFI48; Transcription factor SL1


Immunogen:

Synthesized peptide derived from internal of human TAF1A.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, using TAF1A antibody.