

ADCK2 Antibody



PACO21885

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

ADCK2

Uniprot

Q7Z695

Synonyms:

Uncharacterized aarF domain-containing protein kinase 2; EC 2.7.11; ADCK2; AARF;

Applications:

ELISA, WB

Immunogen:

Synthesized peptide derived from internal of human ADCK2.


Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells and K562 cells, using ADCK2 antibody.