

# ALPK2 Antibody



PACO21888

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

**Protein Background:**

**Gene ID:**

ALPK2

**Uniprot**

Q86TB3

**Synonyms:**

Alpha-protein kinase 2; EC 2.7.11; Heart alpha-protein kinase; ALPK2; HAK

**Applications:**

ELISA, IHC

**Immunogen:**

Synthesized peptide derived from internal of human ALPK2.

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human heart tissue, using ALPK2 antibody.

