

PLXDC1 Antibody



PACO21891

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

PLXDC1

Uniprot

Q8IUUK5

Synonyms:

Plexin domain-containing protein 1; Tumor endothelial marker 7; Tumor endothelial marker 3; TEM3; TEM7

Immunogen:

Synthesized peptide derived from Internal of human PLXDC1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemical analysis of paraffin-embedded human brain tissue using PLXDC1 antibody.

