

TBC1D3 Antibody



PACO21893

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11. Wiemann S. , Genome Res. 11:422-435(2001).

Gene ID:

TBC1D3

Uniprot

Q8IZP1

Synonyms:

TBC1 domain family member 3; Rab GTPase-activating protein PRC17; Prostate cancer gene 17 protein; TRE17 alpha protein; TBC1D3


Immunogen:

Synthesized peptide derived from internal of human TBC1D3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, using TBC1D3 antibody.