

NRIP2 Antibody



PACO21938

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

Protein Background:

Gene ID:

NRIP2

Uniprot

Q9BQI9

Synonyms:

Nuclear receptor-interacting protein 2; NRIP2;

Applications:

ELISA, IHC

Immunogen:

Synthesized peptide derived from N-terminal of human NRIP2.

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using NRIP2 antibody.

