

DCLK3 Antibody



PACO21947

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

DCLK3

Uniprot

Q9C098

Synonyms:

Serine/threonine-protein kinase DCLK3; EC 2.7.11.1; Doublecortin-like and CAM kinase-like 3; Doublecortin-like kinase 3; DCLK3

Immunogen:

Synthesized peptide derived from internal of human DCLK3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using DCLK3 antibody.

