

# ZNF232 Antibody



PACO21986

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

May be involved in transcriptional regulation. Sjoebloom T. , Science 314:268-274(2006).

**Gene ID:**

ZNF232

**Uniprot**

Q9UNY5

**Synonyms:**

Zinc finger protein 232; Zinc finger and SCAN domain-containing protein 11; ZNF232; ZSCAN11;

**Immunogen:**

Synthesized peptide derived from internal of human ZNF232.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ZNF232 antibody.

