

# ZKSCAN5 Antibody



PACO21990

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IF:1:100-1:500

**Protein Background:**

**Gene ID:**

ZKSCAN5

**Uniprot**

Q9Y2L8

**Synonyms:**

Zinc finger protein with KRAB and SCAN domains 5; Zinc finger protein 95 homolog; Zfp-95; ZKSCAN5; KIAA1015

**Immunogen:**

Synthesized peptide derived from internal of human ZNF95.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunofluorescence analysis of A549 cells, using ZNF95 antibody.

