

ADD3 Antibody



PACO22010

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

ADD3

Uniprot

Q9UEY8

Synonyms:

ADDG; ADDL; adducin; gamma; adducin-like protein 70


Immunogen:

Synthesized peptide derived from internal of human ADD3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells and 293 cells, using ADD3 antibody.