

PACO22014

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Belongs to the arrestin family. 2 isoforms of the human protein are produced by alternative splicing.

Gene ID:

ARRDC2

Uniprot

Q8TBH0

Synonyms:

Arrestin domain-containing protein 2; ARRDC2; PP2703;

Immunogen:

Synthesized peptide derived from internal of human ARRDC2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ARRDC2 antibody.

